

## Myasthenia Gravis/Lambert Eaton Syndrome Diagnostic Algorithm

MGLES/83369 Myasthenia Gravis (MG)/Lambert-Eaton Syndrome (LES) Evaluation

Always performed

- P/Q-Type Calcium Channel Antibody, Serum
- N-Type Calcium Channel Antibody, Serum
- Acetylcholine Receptor (Muscle AChR) Binding Antibody, Serum
- Acetylcholine Receptor (Muscle AChR) Modulating Antibodies, Serum
- Striational (Striated Muscle) Antibodies, Serum

If AChR-modulating antibodies:  
loss  $\geq 90\%$   
AND  
Striational antibodies  $\geq 1:60$

Automatically performed\*

- Ganglionic Acetylcholine Receptor (Alpha3) Autoantibody, Serum
- Collapsin Response-Mediator Protein-5-IgG (CRMP-5-IgG) Western Blot, Serum

\* Performed at an additional charge