

**ALDOLASE, S  
#300225**

**EXPLANATION: EXPLANATION:** Effective immediately, Mayo Medical Laboratories will no longer accept specimens drawn in serum gel tubes for test #300225, Aldolase, Serum. Routine validation studies performed by the lab demonstrate that results are optimal when the specimen is collected in a plain, red-top tube.

**SPECIMEN REQUIRED:**

Draw blood in a plain, red-top tube(s) from a fasting patient. **(Serum gel tubes and hemolyzed specimens are not acceptable.)** Spin down and send 1 mL of serum frozen in plastic vial.