

CLOZAPINE, SERUM
#500406

NOTE: This new test will replace 3 tests which will become obsolete on September 16, 2008:

- 500009 – Clozapine, Serum
- 91155 – Clozapine and Metabolite which is referred to NMS Labs
- 91161 – Clozapine and Norclozapine which is referred to MedTox

USEFUL FOR:

- Monitoring patient compliance.
- As an aid in achieving desired therapeutic levels.

METHODOLOGY: Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)

REFERENCE VALUES: Clozapine: >350 ng/mL
Clozapine and Norclozapine: > 450 ng/mL

SPECIMEN REQUIREMENTS: Draw blood in a plain, red-top tube(s). (**Serum gel tube is not acceptable.**) Spin down and send 1.0 mL of serum refrigerated.

CAUTIONS: None

LIST FEE: \$150.30

CPT CODE: 80299 x 2

ANALYTIC TIME: 1 Day

DAY(S) SET UP: Monday through Friday, 3rd shift

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager
Mary Erath, Mayo Medical Laboratories' Technologist Support
Telephone: 800-533-1710

To obtain test set-up information contact MML Client Support or log on to:
MayoMedicalLaboratories.com/test-notifications