

**TEST OBSOLETE**

**NOTIFICATION DATE:** May 7, 2009

**EFFECTIVE DATE:** June 16, 2009

**INDOMETHACIN (INDOCIN), SERUM**  
**#80754**

**EXPLANATION:** Due to low testing volumes, this assay will become obsolete effective June 16, 2009.

**RECOMMENDED ALTERNATIVE TEST:**

**#90525 Miscellaneous National Medical Services (NMS), (NMS test code 2410SP Indomethacin, Serum/Plasma)**

**METHODOLOGY:** High Performance Liquid Chromatography (HPLC)

**SPECIMEN REQUIREMENTS:**

**SUBMIT ONLY 1 OF THE FOLLOWING SPECIMENS:**

**Serum:** Draw blood in a plain, red-top tube(s). (Serum gel tube is not acceptable.) Promptly centrifuge and separate, send 1 mL of serum refrigerate.

**Plasma:** Draw blood in a purple-top (EDTA) tube(s). (Plasma gel tube is not acceptable.) Promptly centrifuge and separate, send 1 mL of EDTA plasma refrigerate.

**Note:** Submission of a serum or plasma gel tubes will result in cancellation.

**Stability:** Refrig 16 days\Ambient 16 days\Frozen 12 months

**LIST FEE:** \$78.20

**CPT CODE:** 87337

**ANALYTIC TIME:** 1 day

**DAY(S) SET-UP:** Monday, Wednesday, Friday

**QUESTIONS:** Contact your Mayo Medical Laboratories' Regional Manager  
Jim Nielsen, Mayo Medical Laboratories' Technologist Support  
Telephone: 800-533-1710