

**CLOZAPINE AND NORCLOZAPINE**  
**#91161**

**EXPLANATION OF CHANGE:** Clozapine and Norclozapine (91161) referred to MedTox labs, will be replaced with new test, Clozapine, Serum (500406).

**UNIT CODE:** 500406

**PUBLISHED NAME:** CLOZAPINE, SERUM

**USEFUL FOR:**

- Monitoring patient compliance.
- As an aid in achieving desired therapeutic levels.

**METHODOLOGY:** Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)

**REFERENCE VALUES:** Clozapine: >350 ng/mL  
Clozapine and Norclozapine: > 450 ng/mL

**SPECIMEN REQUIREMENTS:** Draw blood in a plain, red-top tube(s). (**Serum gel tube is not acceptable.**) Spin down and send 1.0 mL of serum refrigerated.

**CAUTIONS:** None

**LIST FEE:** \$150.30

**CPT CODE:** 80299 x 2

**ANALYTIC TIME:** 1 Day

**DAY(S) SET UP:** Monday through Friday, 3rd shift

**QUESTIONS:** Contact Mayo Medical Laboratories  
Telephone: 800-533-1710

To obtain test set-up information contact MML Client Support or log on to:  
[MayoMedicalLaboratories.com/test-notifications](http://MayoMedicalLaboratories.com/test-notifications)