



MAYO CLINIC
 Mayo Medical Laboratories

TEST DEFINITION

8/27/2008

MML Name/order code cross-reference

Code	Name
91873	Sage IgE

MML MML Test setup information

ORDER CODE	EFF DATE	TC	TITLE	Checking Normals	Print normals (# coded)	Perform Site *
91873	8/25/2008		Sage IgE			
UNIT CODE ALWAYS MESSAGE - [ZEAA]						
TEST PERFORMED BY ESOTERIX ALLERGY AND ASTHMA						
201 SUMMIT VIEW DRIVE						
SUITE 10						
BRENTWOOD, TN 37027						
Transport temp : Refrig\Ambient\Frozen OK						
91873			Sage IgE			
Units: Counts						
24961			Score			
TEST CODE ALWAYS MESSAGE - [ZAIR]						
ZAIR						

Allergen Interpretation Range

Score	IgE Counts	Interpretation
0	0-500	no antibodies detected
0/1	501-750	low equivocal
1	751-1600	equivocal
2	1601-3600	
3	3601-8000	scores 2-6 indicate
4	8001-18000	increasing antibody levels
5	18001-40000	
6	40001+	

*Performing Site Legend

=====
=====

MML Messages used as normals

CODE TEXT

Total of 0 normals codes

*** End of Report ***

To obtain test set-up information contact MML Client Support or log on to:
<http://www.mayomedicallaboratories.com/customer-service/set-up/index.html>