



# TEST DEFINITION

8/19/2008

MML Name/Order Code Cross-Reference  
 Code Name  
 -----  
 91888 Febrile Antibodies Panel

MML MML Test Setup Information

ORDER CODE	EFF DATE	TC	TITLE	Checking Normals	Print normals (# coded)	Perform Site *
-----	-----	-----	-----	-----	-----	-----
91888(profile)			Febrile Antibodies Panel			
90461	2/3/2005		Weil-Felix Test, DA			
			UNIT CODE ALWAYS MESSAGE - [ZFOCUSD]			
			Test Performed by: Focus Diagnostics, Inc.			
			5785 Corporate Avenue			
			Cypress, CA 90630-4750			
			Transport temp : Refrig\Frozen OK\Ambient <4 days OK			
			Serum gel tube stability: Refrig <6 wks\Frozen NO\Ambient NO			
			Worklist : FWEIL			
90461			Proteus Antigen OX19			
			- - - - -			
91024			Proteus Antigen OX2			
			- - - - -			
91025			Proteus Antigen OX K			
			- - - - -			
91889	8/18/2008		Salmonella, Total Antibody, EIA			
			UNIT CODE ALWAYS MESSAGE - [ZFOCUSD]			
			Test Performed by: Focus Diagnostics, Inc.			
			5785 Corporate Avenue			
			Cypress, CA 90630-4750			
			Transport temp : Refrig\Frozen OK\Ambient 7 days OK			
			Worklist : FSALMO			
24955			Salmonella H, Type a			
			- - - - -			
24956			Salmonella H, Type b			
			- - - - -			
24957			Salmonella H, Type d			
			- - - - -			

MML		MML Test setup information			Print normals (# coded)	Perform Site *
ORDER CODE	EFF DATE	TC	TITLE	Checking Normals		
91888	(continued)					
		24958	Salmonella O, Type Vi			
		24959	Salmonella O, Type D			
91847	5/16/2008		Brucella Antibodies (IgG,IgM),Serum			
			UNIT CODE ALWAYS MESSAGE - [ZFOCUSD]			
			Test Performed by: Focus Diagnostics, Inc.			
			5785 Corporate Avenue			
			Cypress, CA 90630-4750			
			Transport temp : Refrig\Frozen OK\Ambient OK			
			Worklist : FBRUA			
		24920	Brucella IgG			
			Units: U			
		24921	Brucella IgM			
			Units: U			

\*Performing Site Legend

MML CODE	TEXT	Messages used as normals

Total of 0 normals codes

\*\*\* End of Report \*\*\*

To obtain test set-up information contact MML Client Support or log on to:  
<http://www.mayomedicallaboratories.com/customer-service/set-up/index.html>