
Reporting Title: CDKN1C Gene, Known Mutation**Performing Location:** Rochester**Specimen Requirements:**

Forms:

1. Molecular Genetics-Congenital Inherited Diseases Patient Information Sheet (Supply T521) in Special Instructions
2. New York Clients-Informed consent is required. Please document on the request form or electronic order that a copy is on file. An Informed Consent for Genetic Testing (Supply T576) is available in Special Instructions.

Specimen must arrive within 96 hours of collection.

Submit only 1 of the following specimens:

Preferred:

Specimen Type: Whole blood

Container/Tube:

Preferred: Lavender top (EDTA) or yellow top (ACD)

Acceptable: Any anticoagulant

Specimen Volume: 3 mL

Collection Instructions:

1. Invert several times to mix blood.

2. Send specimen in original tube.

Specimen Stability Information: Ambient (preferred)/Refrigerated

Specimen Type: Cultured fibroblasts

Container/Tube: T-75 or T-25 flask

Specimen Volume: 1 Full T-75 or 2 full T-25 flasks

Specimen Stability Information: Ambient (preferred)/Refrigerated <24 hours

Due to the complexity of prenatal testing, consultation with the laboratory is required for all prenatal testing. Prenatal specimens can be sent Monday through Thursday and must be received by 5 p.m. CST on Friday in order to be processed appropriately. All prenatal specimens must be accompanied by a maternal blood specimen. Order MCC / Maternal Cell Contamination, Molecular Analysis on both the prenatal and maternal specimens.

Specimen Type: Amniotic fluid

Container/Tube: Amniotic fluid container

Specimen Volume: 20 mL

Specimen Stability Information: Refrigerated (preferred)/Ambient

Specimen Type: Chorionic villi

Container/Tube: 15-mL tube containing 15 mL of transport media

Specimen Volume: 20 mg

Specimen Stability Information: Refrigerated

Acceptable:

Specimen Type: Confluent cultured cells

Container/Tube: T-25 flask

Specimen Volume: 2 Flasks

Collection Instructions: Submit confluent cultured cells from another laboratory.

Specimen Stability Information: Ambient (preferred)/Refrigerated

Specimen Type: Skin biopsy

Container/Tube: Sterile container with any standard cell culture media (eg, minimal essential media, RPMI 1640). The solution should be supplemented with 1% penicillin and streptomycin. Tubes can be supplied upon request (Eagle's minimum essential medium with 1% penicillin and streptomycin [Supply T115]).

Specimen Volume: 4-mm punch

Specimen Stability Information: Refrigerated (preferred)/Ambient

Specimen Type	Temperature	Time
Varies	Varies	

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Type	Reportable
FBC	CG030	Source: • Autopsy • Chorionic Villi • Products of Conception or Stillbirth	Answer List	No
FBC	CG033	Reason for Referral	Plain Text	No
AFC	CG132	Reason For Referral	Plain Text	No

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
33503	Specimen	Alphanumeric		N/A
33504	Specimen ID	Alphanumeric		N/A
33505	Source	Alphanumeric		31208-2
33506	Order Date	Alphanumeric		N/A
33507	Reason for Referral	Alphanumeric		42349-1
33508	Method	Alphanumeric		49549-9
33509	Result	Alphanumeric		N/A
33510	Interpretation	Alphanumeric		N/A
33511	Amendment	Alphanumeric		In Process
33512	Reviewed By	Alphanumeric		N/A

Result ID	Reporting Name	Type	Unit	LOINC®
33513	Release Date	Alphanumeric		N/A

CPT Code: 1 x 81479

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
FBC	Fibroblast Culture for Genetic Test		Profile	No	Yes
AFC	Amniotic Fluid Culture/Genetic Test		Profile	No	Yes
MCC	Maternal Cell Contamination, B	1	81265	No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
FBC	16167	Specimen	Alphanumeric		N/A
FBC	16337	Specimen ID	Alphanumeric		N/A
FBC	CG030	Source	Alphanumeric		31208-2
FBC	16168	Order Date	Alphanumeric		N/A
FBC	CG033	Reason For Referral	Alphanumeric		42349-1
FBC	16172	Method	Alphanumeric		49549-9
FBC	16175	Interpretation	Alphanumeric		69965-2
FBC	16176	Amendment	Alphanumeric		In Process
FBC	16177	Consultant	Alphanumeric		N/A
FBC	16178	Report Date	Alphanumeric		N/A
AFC	16179	Specimen	Alphanumeric		N/A
AFC	16338	Specimen ID	Alphanumeric		N/A
AFC	16793	Source	Alphanumeric		31208-2
AFC	16180	Order Date	Alphanumeric		N/A
AFC	CG132	Reason For Referral	Alphanumeric		42349-1
AFC	16184	Method	Alphanumeric		49549-9

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
AFC	16187	Interpretation	Alphanumeric		69965-2
AFC	16188	Amendment	Alphanumeric		In Process
AFC	16189	Consultant	Alphanumeric		N/A
AFC	16190	Report Date	Alphanumeric		N/A
MCC	23332	Specimen	Alphanumeric		N/A
MCC	23333	Specimen ID	Alphanumeric		N/A
MCC	23334	Source	Alphanumeric		31208-2
MCC	23335	Order Date	Alphanumeric		N/A
MCC	23338	Result	Alphanumeric		40704-9
MCC	23337	Amendment	Alphanumeric		40704-9
MCC	23341	Reviewed By:	Alphanumeric		N/A
MCC	23339	Release Date	Alphanumeric		N/A

Reference Values:

An interpretive report will be provided.