Creatinine, Urine



**Reporting Title:** Creatinine, U **Performing Location:** Rochester

### **Specimen Requirements:**

Container/Tube: Plastic, 6-mL tube (Supply T465)

Specimen Volume: 5 mL Collection Instructions:

- 1. Collect urine for 24 hours.
- 2. No preservative. Additional Information:

Thymol Yes

- 1. 24 Hour volume is required.
- 2. See Urine Preservatives in Special Instructions for multiple collections.

Urine Preservative Collection Options
Ambient Preferred
Refrigerated Yes
Frozen Yes
6N HCI Yes
50% Acetic Acid Yes
Na2CO3 Yes
Toluene Yes
6N HNO3 Yes
Boric Acid Yes

Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	30 days
	Ambient	30 days
	Frozen	30 days

## Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
СТИ	TM10	Collection Duration	Plain Text	Yes
СТИ	VL8	Urine Volume	Plain Text	Yes



# **Test Definition: CTU**

Creatinine, Urine

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
CRE_A	Creatinine, U	Numeric	mg/spec	2162-6
TM10	Collection Duration	Numeric	h	13362-9
VL8	Urine Volume	Numeric	mL	3167-4
NCTU_	Creatinine Conc	Numeric	mg/dL	2161-8

**CPT Code:** 1 × 82570

#### **Reference Values:**

15-25 mg/kg of body weight/24 hours Reported in unit of mg/specimen