

## **Test Definition: FASS**

Antidepressant Drug Screen, S/P

Reporting Title: Antidepressant Drug Screen, S/P

Performing Location: Medtox Laboratories, Inc

## **Specimen Requirements:**

Submit only 1 of the following:

#### Serum

Draw blood in a plain, red-top tube(s). (Serum gel tube is not acceptable.) Spin down and send 3 mL of serum refrigerate in plastic vial

#### Plasma

Draw blood in a green-top (sodium heparin) tube(s). (Plasma gel tube is not acceptable.) Spin down and send 3 mL of plasma refrigerate in plastic vial.

Specimen Type	Temperature	Time	
Varies	Refrigerated (preferred)		
	Ambient	72 hours	
	Frozen		

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z1115	Submitted Specimen Type:	Alphanumeric		In Process
Z1116	Antidepressant Screen	Alphanumeric		In Process
Z1117	Antidepressant Screen Confirmation	Alphanumeric		In Process

## **Components:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
Billing only	Column chromatography	1	82541		
Billing only	Antidepressant Drug Screen, S/P	1	80101		



# **Test Definition: FASS**

Antidepressant Drug Screen, S/P

#### **Reference Values:**

Antidepressant screen includes the analysis for:
Amytriptyline, Clomipramine and Desmethylclomipramine,
Cyclobenzaprine, Desipramine, Doxepin and Desmethyldoxepin,
Fluoxetine and Norfluoxetine, Imipramine, Maprotiline,
Nortriptyline, Paroxetine, Protriptyline, Sertraline and
Desmethylsertraline, Trimipramine.

Test Performed By: Medtox Laboratories

402 W. County Road D St. Paul, MN 55112