

Test Definition: MEX

Mexiletine, Serum

Reporting Title: Mexiletine, S Performing Location: Rochester

Specimen Requirements:

Container/Tube: Red top Specimen Volume: 3 mL

Collection Instructions: Draw blood after patient has been receiving mexiletine for at least 3 days and just before

administration of the next dose.

Specimen Type	Temperature	Time
Serum Red	Refrigerated (preferred)	14 days
	Frozen	14 days
	Ambient	7 days

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
9245	Mexiletine, S	Numeric	mcg/mL	3819-0

CPT Code: 1 × 80299

Reference Values:

Therapeutic concentration: 0.8-2.0 mcg/mL (trough value)

Toxic concentration: >2.0 mcg/mL (trough value)