

## REPORTING CHANGE

**NOTIFICATION DATE:** June 29, 2010 **EFFECTIVE DATE:** July 27, 2010

## SYNTHETIC GLUCOCORTICOID SCREEN, SERUM #81031

**EXPLANATION:** The cutoff concentrations of these analytes will change due to the precision of the assay. The concentration values will change from ug/dL to mcg/dL.

**NOTE**: Beclomethasone Dipropionate is being removed from this assay.

Analyte Name	Current Cutoff Concentration	New Cutoff Concentration
Betamethasone	0.05 ug/dL	0.10 mcg/dL
Budesonide	0.04 ug/dL	$0.10~\mathrm{mcg/dL}$
Dexamethasone	0.06 ug/dL	0.10 mcg/dL
Fludrocortisone	0.06 ug/dL	0.10 mcg/dL
Flunisolide	0.03 ug/dL	0.10 mcg/dL
Fluorometholone	0.04 ug/dL	0.10 mcg/dL
Fluticasone propionate	0.04 ug/dL	0.10 mcg/dL
Megestrol acetate	0.03 ug/dL	0.10 mcg/dL
Methylprednisolone	0.06 ug/dL	0.10 mcg/dL
Prednisolone	0.05 ug/dL	0.10 mcg/dL
Prednisone	0.05 ug/dL	0.10 mcg/dL
Triamcinolone acetonide	0.03 ug/dL	0.10 mcg/dL

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or Jeremy Caddell, MML Laboratory Technologist Resource Coordinator Telephone: 800-533-1710