

REPORTING AND METHODOLOGY CHANGE

NOTIFICATION DATE: August 16, 2010 EFFECTIVE DATE: September 8, 2010

Mannose-Binding Lectin, Serum #81051

EXPLANATION: Beginning September 8th, this test will be converted from a manual assay to an automated testing platform. As a result, the reference values and reportable range for this assay will change. All other characteristics of this test will remain the same.

CURRENT REFERENCE VALUES: > or = 5.0 ng/mL

NEW REFERENCE VALUES: > or = 7.8 ng/mL

CURRENT REPORTING RANGE: 5.0 - 4000 ng/mL

NEW REPORTING RANGE: 5.0 - 400 ng/mL

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager Brad Bellin, Mayo Medical Laboratories' Technical Support Telephone: 800-533-1710