

REFERENCE VALUE CHANGE

NOTIFICATION DATE: September 28, 2010 **EFFECTIVE DATE:** October 27, 2010

MANGANESE, RANDOM, URINE #86154

EXPLANATION: Reference values will be added to Manganese, Random, Urine, test #86154 based on a Mayo Medical Laboratories Metals Lab reference value study.

CURRENT REFERENCE VALUE: No established reference value

NEW REFERENCE VALUE:

0-15 years: Not Established > or = 16 years: 0.1-0.9 mcg/L

NOTE: High concentrations of gadolinium and iodine are known to interfere with most metals tests. If either gadolinium- or iodine-containing contrast media has been administered, a specimen cannot be collected for 96 hours.

NOTE: No changes in CPT or fee will occur.

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or Jim Nielsen, MML Laboratory Technologist Resource Coordinator Telephone: 800-533-1710