## MAYO CLINIC Mayo Medical Laboratories

## **TEST DEFINITION**

3/17/2011

ORDER EFF PERFORM CODE DATE TC TITLE PRINT NORMALS (# CODED) SITE \* CHECKING NORMALS 200814 3/11/2011 Breast Carcinoma Assoc Ag(CA 27.29) NEL UNIT CODE ALWAYS MESSAGE - [ZCA2729A] THE TESTING METHOD IS A CHEMILUMINOMETRIC IMMUNOASSAY MANUFACTURED BY SIEMENS AND PERFORMED ON THE SIEMEN'S ADVIA CENTAUR. VALUES OBTAINED WITH DIFFERENT ASSAY METHODS OR KITS MAY BE DIFFERENT AND CANNOT BE USED INTERCHANGEABLY. TEST RESULTS CANNOT BE INTERPRETED AS ABSOLUTE EVIDENCE FOR THE PRESENCE OR ABSENCE OF MALIGNANT DISEASE. TRANSPORT TEMP: FROZEN\REFRIG < 48 HOURS OK\AMBIENT NO 200814 Breast Carcinoma Assoc Ag(CA 27.29) UNITS: U/ML NO SEX NO AGE ; #N81413N 1D TO 17Y : ; #N81413Y18 18Y TO 150Y : 0.0-38.0 i <= 38.0MALE No age : ; #N81413N 1D TO 17Y : ; #N81413Y18 18Y TO 150Y : 0.0-38.0 ; <=38.0 FEMALE No age : ; #N81413N 1D to 17Y : ; #N81413Y18 18Y TO 150Y : 0.0-38.0 ; <=38.0 \*Performing Site Legend NEL Mayo Medical Laboratories New England LAB DIRECTOR: LYNN A. CHERYK, PH.D. 160 DASCOMB ROAD ANDOVER, MA 01810 \_\_\_\_\_\_

CODE TEXT

\*\*\* END OF REPORT \*\*\*

\_\_\_\_\_

N81413N REFERENCE VALUES HAVE NOT BEEN

ESTABLISHED FOR PATIENTS THAT ARE LESS

THAN 18 YEARS OF AGE.

<=38.0 (>= 18<sub>Y</sub>)

-----

N81413Y18 REFERENCE VALUES HAVE NOT BEEN

ESTABLISHED FOR PATIENTS THAT ARE LESS

THAN 18 YEARS OF AGE.

-----

Total of 2 normals codes

\*\*\* END OF REPORT \*\*\*