

SPECIMEN REMINDER

NOTIFICATION DATE: September 22, 2011

EFFECTIVE DATE: Immediately

PRENATAL PROFILE III
#31794

EXPLANATION: This custom profile test requires both EDTA whole blood and serum. Effective immediately, we will no longer be sharing whole blood EDTA specimens for the CBC and blood bank testing. Consequently, please send at least two EDTA whole blood tubes when ordering this profile. In cases where only one EDTA whole blood tube is received, the CBC test (#9109, CBC with Differential, Blood) will be canceled.

NOTE: This change may impact test set-up information and could require a change to file definition. Please review the Test Set-Up information for specifics.

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or
Richard R. Einerson, Mayo Medical Laboratories' Technical Support
Telephone: 800-533-1710