



MAYO CLINIC
Mayo Medical Laboratories

1-800-533-1710

TEST UPDATE

REFERRAL

NOTIFICATION DATE: April 4, 2012

EFFECTIVE DATE: Immediately

LATEX SPECIFIC IGE (RIA) PANEL

Test ID: FLTX

Secondary ID: 57118

EXPLANATION OF CHANGE: Test ID FLTX, referred to Viracor-IBT Laboratories, will reflect the following update effective immediately.

CURRENT SPECIMEN REQUIREMENTS: Draw blood in a plain, red-top tube(s). Spin down and send 0.5 mL of serum.

NEW SPECIMEN REQUIREMENTS: Draw blood in a plain, red-top tube(s). Spin down and send 2 mL of serum refrigerate in a plastic vial.

NOTE: This change may impact test set-up information and could require a change to file definition. Please review the Test Set-Up information for specifics at <http://www.mayomedicallaboratories.com/test-notifications/index.html>.

QUESTIONS: Contact Mary Erath, MML Laboratory Technologist Resource Coordinator
Telephone: 800-533-1710