TRANSPORT TEMPERATURE



CHANGE REFERRAL

NOTIFICATION DATE: April 18, 2012 **EFFECTIVE DATE:** Immediately

MVISTA HISTOPLASMA ANTIGEN EIA

Test ID: FHST Secondary ID: 91957

EXPLANATION: Test ID FHST, referred to MiraVista Diagnostics, will reflect the following update effective immediately.

CURRENT TRANSPORT TEMPERATURE:

Specimen Type	Temperature	Time
Varies	Refrigerated (preferred)	30 days
	Ambient	5 days
	Frozen	

NEW TRANSPORT TEMPERATURE:

Specimen Type	Temperature	Time
Varies	Refrigerated (preferred)	14 days
	Ambient	14 days
	Frozen	

QUESTIONS: Contact Mary Erath, MML Laboratory Technologist Resource Coordinator Telephone: 800-533-1710