

## **Test Definition: FAAST**

Adenovirus Antigen, EIA Stool

Reporting Title: Adenovirus Antigen, EIA (Stool)

Performing Location: Focus Diagnostics,

## **Specimen Requirements:**

Submit only 1 of the following specimens:

Stoo

Collect 2 grams or 2 mL fresh unpreserved stool in sterile container. Ship frozen.

Rectal Swab

1 rectal swab. Ship frozen.

| Specimen Type | Temperature        | Time     |
|---------------|--------------------|----------|
| Varies        | Frozen (preferred) | 14 days  |
|               | Refrigerated       | 72 hours |

## **Result Codes:**

| Result ID | Reporting Name          | Туре         | Unit | LOINC® |
|-----------|-------------------------|--------------|------|--------|
| Z2275     | Adenovirus Antigen, EIA | Alphanumeric |      | 5825-5 |

**CPT Code:** 1 × 87301

## **Reference Values:**

Not Detected

A positive test result in a patient with diarrhea strongly suggests that Adenovirus is the cause of the gastroenteritis. However, Adenovirus may be shed in stool asymptomatically as long as 18 months post-infection. Enterocolitis due to Staphylococcus aureus infection may cause a false positive result with this test due to high levels of bacterial protein A.

Test Performed by: Focus Diagnostics, Inc 5785 Corporate Ave. Cypress, CA 90630