
Reporting Title: Enoxaparin Dependent Platelet Ab
Performing Location: Bloodcenter Of Wisc

Specimen Requirements:

Draw blood in a plain, red-top tube. (Serum gel tube is acceptable.) Spin down and remove serum from clot. Ship 5 mL of serum refrigerated in a plastic vial.

Note: Date of birth required.

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	7 days
	Frozen	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z0126	enoxaparin Low Dose	Alphanumeric	%	
Z0127	enoxaparin High Dose	Alphanumeric	%	
Z0128	enoxaparin Result	Alphanumeric		

CPT Code: 1 x 86022

Reference Values:

Negative

Test Performed by: BloodCenter of Wisconsin
638 N. 18th Street
Milwaukee, WI 53233-2121