

# **Test Definition: BRUGM**

Brucella Antibody Screen, IgG & IgM, Serum

Reporting Title: Brucella Ab Screen, IgG and IgM, S

**Performing Location:** Rochester

## **Specimen Requirements:**

Container/Tube: Preferred: Red top Acceptable: Serum gel Specimen Volume: 2 mL

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Frozen	14 days

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
23685	Brucella Ab Screen, IgG, S	Alphanumeric		
23686	Brucella Ab Screen, IgM, S	Alphanumeric		
23687	Interpretation	Alphanumeric		

### **Components:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
Billing only	Brucella Ab Screen, IgM, S	1	86622		
Billing only	Brucella Ab Screen, IgG, S	1	86622		

### **Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
BRUTA	Brucella Ab, Agglutination, S	1	86622	No	Yes



# **Test Definition: BRUGM**

Brucella Antibody Screen, IgG & IgM, Serum

### **Result Codes for Reflex Tests:**

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
BRUTA	8112	Brucella Ab, Agglutination, S	Alphanumeric		40614-0

#### **Reference Values:**

IgG SCREEN

Negative (reported as positive, negative, or equivocal)

IgM SCREEN

Negative (reported as positive, negative, or equivocal)