

## **TEST DOWN NOTIFICATION**

MML NEW ENGLAND

**NOTIFICATION DATE:** August 6, 2012 **EFFECTIVE DATE:** Immediately

## **INSULIN-LIKE GROWTH FACTOR 1, SERUM**

Test ID: IGF1 Secondary ID: 200009

**EXPLANATION:** Due to unavailability of reagent from the manufacturer, this test will become non-orderable effective immediately.

**NOTE:** Once this has been resolved, the test code may be reactivated.

RECOMMENDED ALTERNATE TEST: Test ID: FILGF

Test Name: Insulin-Like Growth Factor 1 (IGF-1)

Referred to: Esoterix Endocrinology

**METHODOLOGY:** Blocking RIA after acid alcohol extraction

## **REFERENCE VALUES:**

	Term		Preterm	
Age	Range	Mean	Range	Mean
Birth	15-109	59	21-93	51
2m	15-109	55	23-163	81
4m	7-124	50	23-171	74
6m	7-93	41	15-132	61
12m	15-101	56	15-179	77
	Male		Female	
Age	Range	Mean	Range	Mean
1-2y	30-122	76	56-144	100

Pubertal ranges are available in the Endocrinology Syllabus at Esoterix.com

	Male		Female	
Age	Range	Mean	Range	Mean
19-20y	281-510	371	217-475	323
21-30y	155-432	289	87-368	237
31-40y	132-333	226	106-368	225
41-50y	121-237	160	118-298	205
51-60y	68-245	153	53-287	172
61-70y	60-220	132	75-263	180
71-80y	36-215	131	54-205	156

**SPECIMEN REQUIREMENTS:** 

**Container/Tube:** Serum Red Top

**Collection Instructions:** Spin down and separate within one hour of collection, freeze

immediately. Ship 0.5 mL frozen in a plastic vial.

## **SPECIMEN STABILITY INFORMATION:**

Specimen Type	Temperature	Time
Serum	Frozen (preferred)	
	Ambient	48 hours
	Refrigerated	48 hours

**FEE:** \$ 95.00 **CPT CODE**: 84305

**DAY(S) SET UP:** Monday - Friday **ANALYTIC TIME:** 3 – 5 days

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or Rita Baird, MML Laboratory Technologist Resource Coordinator Telephone: 800-533-1710