

PROCAINAMIDE, SERUM

Test ID: PA

Secondary ID: 8683

NOTE: Effective with this change, all orders placed for this test must be submitted in the alpha-numeric character Test ID format instead of the numeric Secondary ID. Result codes may also be subject to change. Please review Test Set-Up and Conversion Mapping information at: <http://www.mayomedicallaboratories.com/test-notifications/index.html>

EXPLANATION: The reference values for Procainamide, Serum are changing.

CURRENT REFERENCE VALUES:

Test	Reference Values	Toxic Concentration	Units
Procainamide	4.0-8.0 (Therapeutic concentration)	>=16.0	mcg/mL
N-Acetyl Procainamide	4.0-8.0 (Therapeutic concentration)	>=16.0	mcg/mL
Procainamide and N-Acetylprocainamide	<=30.0 (total of both) (Therapeutic concentration)	>=30.0	mcg/mL

NEW REFERENCE VALUES:

Test	Reference Values	Toxic Concentration	Units
Procainamide	4.0-8.0 (Therapeutic concentration)	>10.0	mcg/mL
N-Acetyl Procainamide	10.0-20.0 (Therapeutic concentration)	>40.0	mcg/mL
Procainamide + NAPA			mcg/mL

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or Jim Nielsen, MML Laboratory Technologist Resource Coordinator
 Telephone: 800-533-1710