

---

**Reporting Title:** PSA, Ultrasensitive  
**Performing Location:** Labcorp Burlington

**Specimen Requirements:**

Draw blood in a red-top tube or serum gel tube. Spin down and send 0.8 mL of serum in plastic screw-capped tube, shipped frozen.

Specimen Type	Temperature	Time
Serum	Frozen (preferred)	14 days
	Ambient	7 days
	Refrigerated	7 days

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
Z0983	PSA, Ultrasensitive	Numeric	ng/mL	35741-8

**CPT Code:** 1 x 84153

**Reference Values:**

0.0 00 - 4.000 ng/mL

Roche ECLIA methodology

According to the American Urological Association, Serum PSA should decrease and remain at undetectable levels after radical prostatectomy. The AUA defines biochemical recurrence as an initial PSA value 0.200 ng/mL or greater followed by a subsequent confirmatory PSA value 0.200 ng/mL or greater.

Values obtained with different assay methods or kits cannot be used interchangeably. Results cannot be interpreted as absolute evidence of the presence or absence of malignant disease.

Test Performed by: LabCorp Burlington  
1447 York Court  
Burlington, NC 27215