

## **Test Definition: EXHB**

Hematologic Disorders, DNA/RNA Extract and Hold, Blood

Reporting Title: Heme Extract and Hold, B

**Performing Location: Rochester** 

## **Specimen Requirements:**

Specimen must arrive within 72 hours of draw.

Collection Container/Tube: Preferred: Lavender top (EDTA)

Acceptable: ACD
Specimen Volume: 4 mL
Collection Instructions:

Invert several times to mix blood.
 Send specimen in original tube.

Additional Information: Date of draw is required.

Specimen Type	Temperature	Time	
Whole blood	Ambient (preferred)	72 hours	
	Refrigerated	72 hours	

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
60562	DNA and RNA Extraction	Alphanumeric		In Process