

Reporting Title: OVA1**Performing Location:** Quest Diag-chantill**Specimen Requirements:**

Draw blood in a plain red-top tube(s). (Serum gel tube is acceptable.) Spin down and submit 2.2 mL of serum frozen in a plastic vial.

Note: OVA1 Test is not to be used if patient is known to have rheumatoid factor of > or = 250 IU/L. Triglycerides above 450 mg/dL may interfere with the assay and should not be used. The test is not indicated for patients with diagnosis of malignancy within the last 5 years.

Specimen Type	Temperature	Time
Serum	Frozen (preferred)	63 days
	Refrigerated	5 days

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FOVA1	OVA1	Numeric		In Process

CPT Code Information:

84999

Reference Values:

OVA1 Test Value

Interpretive Criteria for Probability of Malignancy

Low High

Premenopausal <5.0 > = 5.0

Postmenopausal <4.4 > = 4.4

The OVA1 Test is based on a proprietary FDA-cleared software device that combines into a single numerical result; the results of Beta-2 microglobulin, Transferrin, Transthyretin, Apolipoprotein A1 and CA 125 II (Roche Elecsys). These tests should be ordered independently if clinically indicated. The OVA1 Test is an aid to further assess the likelihood that malignancy is present when the physician's independent clinical and radiological evaluation does not indicate malignancy.