



MAYO
Mayo Medical Laboratories
1-800-533-1710

NEW TEST ANNOUNCEMENT
REFERRAL

NOTIFICATION DATE: November 14, 2012
EFFECTIVE DATE: December 10, 2012

FLT3 MUTATION TESTING
Test ID: FFLT3

EXPLANATION: The following new orderable, referred to Laboratory for Personalize Molecular Medicine, will be available December 10, 2012.

REFERRAL LAB CODE: FLT3

METHODOLOGY: Multiplex PCR

REFERENCE VALUES: Negative

SPECIMEN REQUIREMENTS:

Submit only 1 of the following specimens:

- Note:
1. If also ordering a Mayo test, send separate specimens
 2. Clinical reason for referral (AML or CML) is required
 3. Specimen source is required
 4. Date of collection is required

Blood

Draw blood in a lavender-top (EDTA) or yellow-top (ACD) tube(s). Send 4 mL EDTA or ACD whole blood refrigerate.

Bone Marrow

Collect bone marrow in a lavender-top (EDTA) or yellow-top (ACD) tube(s). Send 3 mL EDTA or ACD bone marrow refrigerate.

DNA

1 ug previously isolated DNA shipped refrigerate.

SPECIMEN STABILITY INFORMATION:

Specimen Type	Temperature	Time
Varies	Refrigerated (preferred)	7 days
	Ambient	72 hours

FEE: \$325.00

CPT Code: 83891, 83898 x2, 83892, 83909 x2, 83912

DAY(S) SET UP: Varies

ANALYTIC TIME: 3 – 5 days

QUESTIONS: Contact Mary Erath, MML Laboratory Technologist Resource Coordinator
 Telephone: 800-533-1710