

**Reporting Title: DRVVT**

**Performing Location: Esoterix Coagulatio**

**Specimen Requirements:**

Draw blood in a light blue-top (Sodium citrate) tube(s). Spin down and send 1 mL citrated plasma frozen in a plastic vial.

Note: Separate specimens must be submitted when multiple tests are ordered.

Specimen Type	Temperature	Time
Plasma Na Cit	Frozen (preferred)	30 days
	Refrigerated	4 hours

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
Z0214	DRVVT Screen Seconds	Numeric	sec	6303-2
Z0215	DRVVT Confirm Seconds	Alphanumeric	sec	57838-5
Z0216	DRVVT Ratio	Numeric	ratio	15359-3

**Components:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
FSCRN	DRVVT Screen Seconds	1	85613	Yes	No
FCONF	DRVVT Confirm Seconds	1	85613	No	No
FRATI	DRVVT Ratio			No	No

**Reference Values:**

DRVVT Screen Seconds  
<55.1

DRVVT Confirm Seconds  
No reference range provided.

DRVVT Ratio  
<1.3

DRVVT Ratio may be factitiously elevated when patient is on anti-vitamin K therapy. Esoterix Coagulation recommends performing additional confirmatory tests prior to diagnosing or repeating DRVVT when patient is off anti-vitamin K therapy.

Test Performed By: Esoterix Coagulation  
8490 Upland Dr.  
Suite 100  
Englewood, CO 80112