

TEST DOWN NOTIFICATION

NOTIFICATION DATE: November 30, 2012

EFFECTIVE DATE: Immediately

BAKER'S YEAST, IGE

Test ID: BYST Secondary ID: 82759

EXPLANATION: Due to unavailability of reagents, BYST, Baker's Yeast, IgE assay will be down until further notice effective immediately. Orders for this test received 12/4/2012 and thereafter, including profiles containing this test, will be changed to FYSTB Brewers Yeast (Saccaromyces spp) IgE, referred to Viracor-IBT Laboratories. Clients with pending orders for this testing will be contacted and offered two options:

- Forward test to alternate lab
- Cancel test

RECOMMENDED ALTERNATE TEST: FYSTB Brewers Yeast (Saccaromyces spp) IgE, referred to Viracor-IBT Laboratories.

METHODOLOGY: Radioimmunoassay (RIA)

REFERENCE VALUES:

Units kU/L Reference Range <0.35

SPECIMEN REQUIREMENTS: Draw blood in a plain, red-top tube(s) or a serum gel tube(s). Spin down and send 0.5 mL of serum refrigerated.

Specimen Stability Information

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	28 days
	Ambient	28 days
	Frozen	365 days

LIST FEE: \$34.36 **CPT CODE:** 86003

DAY(S) SET UP: Monday through Friday **ANALYTIC TIME:** 2-3 days

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or Greg Renkly, MML Laboratory Technologist Resource Coordinator Telephone: 800-533-1710