

Test Definition: HEXU

loxehol, Urine

Reporting Title: Iohexol, U **Performing Location:** Rochester

Specimen Requirements:

Container/Tube: Plastic, 6-mL tube (Supply T465)

Specimen Volume: 5 mL Collection Instructions:

1. Collect a timed urine specimen.

2. Patient should be fasting for 4 hours; 2 hours for diabetic patients.

Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	7 days
	Frozen	35 days

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
61712	Iohexol, U	Numeric	mcg/mL	In Process

CPT Code Information:

83789-LC/MS/MS