

**SPECIMEN  
STABILITY/ACCEPTABILITY  
CHANGE****MML JACKSONVILLE****NOTIFICATION DATE:** December 13, 2012**EFFECTIVE DATE:** December 18, 2012**T- AND B-CELL QUANTITATION BY FLOW CYTOMETRY**  
Test ID: 800295

**EXPLANATION:** Effective December 18, 2012, Refrigerated specimens will no longer be acceptable for the T and B Cell Quantitation by Flow Cytometry #800295.

**CURRENT SPECIMEN STABILITY INFORMATION:**

Specimen Type	Temperature	Time
Whole blood	Ambient (preferred)	72 hours
	Refrigerated	72 hours

**NEW SPECIMEN STABILITY INFORMATION:**

Specimen Type	Temperature	Time
Whole blood	Ambient	72 hours

**QUESTIONS:** Contact your Mayo Medical Laboratories' Regional Manager or  
Greg Renkly, MML Laboratory Technologist Resource Coordinator  
Telephone: 800-533-1710