

**METHOD/FILE DEFINITION
CHANGE**

NOTIFICATION DATE: February 19, 2013

EFFECTIVE DATE: March 20, 2013

HLA-B 5701 GENOTYPE, ABACAVIR HYPERSENSITIVITY, SALIVA
Test ID: HL57O
Secondary ID: 60347

NOTE: Effective with this change, all orders placed for this test must be submitted in the alpha-numeric character Test ID format instead of the numeric Secondary ID. Result codes may also be subject to change. Please review Test Set-Up and Conversion Mapping information at: <http://www.mayomedicallaboratories.com/test-notifications/index.html>

EXPLANATION: Effective March 20, 2013, Test HL57O will be performed using a different methodology.

CURRENT METHODOLOGY: Allele-Specific Polymerase Chain Reaction (PCR) Amplification

NEW METHODOLOGY: Qualitative allele-specific real-time Polymerase Chain Reaction (PCR)

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or
Shirley Pokorski, MML Laboratory Technologist Resource Coordinator
Telephone: 800-533-1710