

Reporting Title: Anti-IgA**Performing Location: VIRACOR-IBT LABORATORIES, INTERFACE****Specimen Requirements:**

Draw blood in a plain, red-top tube(s) or serum gel tube(s). Spin down and send 1 mL of serum refrigerate in a plastic vial.

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Frozen	365 days
	Ambient	48 hours

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FIGA	Anti-IgA	Alphanumeric		In Process

CPT Code: 1 x 83520**Reference Values:**

<113 U/mL

Patients with IgG antibodies against IgA may suffer from anaphylactoid reactions when given IVIG that contains small quantities of IgA. In one study (Clinical Immunology 2007; 122:156) five out of eight patients with IgG anti-IgA antibodies developed anaphylactoid reactions when IVIG was administered.

Test Performed by: Viracor-IBT Laboratories
1001 NW Technology Dr
Lee's Summit, MO 64086