

Test Definition: FGLUT

Gluten IgG

Reporting Title: Gluten IgG

Performing Location: VIRACOR-IBT LABORATORIES, INTERFACE

Specimen Requirements:

Draw blood in a plain, red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	7 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
FGLUT	Gluten IgG	Alphanumeric		In Process

CPT Code: 1 × 86001

Reference Values:

< 2 mcg/mL

The reference range listed on the report is the lower limit of quantitation for the assay. The clinical utility of food-specific IgG tests has not been established. These tests can be used in special clinical situations to select foods for evaluation by diet elimination and challenge in patients who have food-related complaints. It should be recognized that the presence of food-specific IgG alone cannot be taken as evidence of food allergy and only indicates immunologic sensitization by the food allergen in question. This test should only be ordered by physicians who recognize the limitations of the test.

Test Performed by: Viracor-IBT Laboratories 1001 NW Technology Dr. Lee's Summit, MO 64086