

## **TEST OBSOLETE**

**NOTIFICATION DATE:** March 26, 2013 **EFFECTIVE DATE:** April 27, 2013

## HLA CLASS II MOLECULAR TYPING DISEASE ASSOCIATION

Test ID: DISII Secondary ID: 32864

**EXPLANATION**: Effective April 27, 2013, due to changes in HLA nomenclature, MML will obsolete test DISII, HLA Class II Molecular Typing Disease Association. Alternatively, clients may order code DIS2, HLA Class II Molecular Typing Disease Association, which incorporates a new reporting format with nomenclature adopted by the WHO Nomenclature Committee for Factors of the HLA System.

**RECOMMENDED ALTERNATIVE TEST:** DIS2, HLA Class II Molecular Typing Disease Association

**USEFUL FOR:** Determining HLA Class II compatibility on specimens for bone marrow and solid organ transplant candidates and their donors.

**REFERENCE VALUES:** Not applicable

**SPECIMEN REQUIREMENTS:** 

**Container/Tube:** Yellow top (ACD solution B)

**Specimen Volume:** 6 mL

**Collection Instructions:** Do not transfer blood to other containers.

Additional Information: Specimen acceptability is based on extracted DNA concentration and

not sample age.

## **CPT CODES:**

81375-HLA Class II typing, low resolution (e.g., antigen equivalents); HLA-DRB1/3/4/5 and -DQB1

**ANALYTIC TIME:** 5-9 days **DAY(S) SET-UP:** Monday through Friday; Monday

through Friday; 7:30 a.m.-5 p.m.

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or Richard R. Einerson, MML Laboratory Technologist Resource Coordinator Telephone: 800-533-1710