## **Test Definition: JAKXM**

JAK2 Exon 12 and Other Non-V617F Mutation Detection, Bone Marrow

Reporting Title: JAK2 Exon 12 Mutation Detection, BM

Performing Location: Rochester

## **Specimen Requirements:**

Specimen must arrive within 72 hours of collection.

Container/Tube:

Preferred: Lavender top (EDTA)

Acceptable: ACD Specimen Volume: 2 mL Collection Instructions:

1. Invert several times to mix bone marrow.

2. Send specimen in original tube.

Additional Information: Collection date is required.

Forms: Hematopathology Patient Information Sheet (Supply T676) in Special Instructions

Specimen Type	Temperature	Time
Bone Marrow	Ambient (preferred)	72 hours
	Refrigerated	72 hours

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
20250	Final Diagnosis:	Alphanumeric		34574-4

**CPT Code:** 1 × 81403

## **Reference Values:**

An interpretive report will be provided.