

# **Test Definition: FECVA**

**Echovirus Antibodies** 

Reporting Title: Echovirus Antibodies

Performing Location: ARUP Laboratories, Interface

## **Specimen Requirements:**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Separate from cells within 2 hours of collection. Send 1 mL of serum refrigerate in a plastic vial.

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	7 days
	Frozen	365 days
	Ambient	48 hours

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z3287	Echovirus Antibody Type 6	Alphanumeric		In Process
Z3288	Echovirus Antibody Type 7	Alphanumeric		In Process
Z3289	Echovirus Antibody Type 9	Alphanumeric		In Process
Z3290	Echovirus Antibody Type 11	Alphanumeric		In Process
Z3291	Echovirus Antibody Type 30	Alphanumeric		In Process

### **CPT Code Information:**

86658 x 5

#### **Reference Values:**

	Ref Interval
Echovirus Antibody Type 6	<1:10
Echovirus Antibody Type 7	<1:10
Echovirus Antibody Type 9	<1:10
Echovirus Antibody Type 11	<1:10
Echovirus Antibody Type 30	<1:10

Interpretive Information: Echovirus Antibody Type 30:

Single positive antibody titers of greater than or equal to 1:80 may indicate past or current infection. Sero-conversion or an increase in titers between acute and convalescent sera of at least fourfold is considered strong evidence of current or recent infection.



# **Test Definition: FECVA**

**Echovirus Antibodies** 

Test Performed by: ARUP Laboratories

500 Chipeta Way

Salt Lake City, UT 84108