

Test Definition: FBEFG

Beef IgG

Reporting Title: Beef IgG

Performing Location: VIRACOR-IBT LABORATORIES, INTERFACE

Specimen Requirements:

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature Time	
Serum	Refrigerated (preferred)	28 days
	Frozen	365 days
	Ambient	7 days

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
FBEFG	Beef IgG	Alphanumeric		In Process

CPT Code Information:

86001

Reference Values:

<2 mcg/mL

The reference range listed on the report is the lower limit of quantitation for the assay. The clinical utility of food-specific IgG tests has not been established. These tests can be used in special clinical situations to select foods for evaluation by diet elimination and challenge in patients who have food-related complaints. It should be recognized that the presence of food-specific IgG alone cannot be taken as evidence of food allergy and only indicated immunologic sensitization by the food allergen in question. This test should only be ordered by physicians who recognize the limitations of the test.

Test Performed by: Viracor-IBT Laboratories

1001 NW Technology Dr. Lee's Summit, MO 64086