

SPECIMEN STABILITY CHANGE

NOTIFICATION DATE: May 1, 2013

EFFECTIVE DATE: Immediately

Lactate Dehydrogenase (LD), Body Fluid

Test ID: LDBF

Secondary ID: 8022

EXPLANATION: Effective immediately, frozen specimens for LDBF, Lactate Dehydrogenase (LD), Body Fluid, are no longer acceptable and will be rejected.

OLD SPECIMEN STABILITY INFORMATION:

Specimen Type	Temperature	Time
Body Fluid	Ambient (preferred)	7 days
	Frozen	30 days
	Refrigerated	48 hours

NEW SPECIMEN STABILITY INFORMATION:

Specimen Type	Temperature	Time
Body Fluid	Ambient (preferred)	7 days
	Refrigerated	48 hours

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or
Richard R. Einerson, MML Laboratory Technologist Resource Coordinator
Telephone: 800-533-1710