

RUBEOLA (MEASLES) ANTIBODIES SCREEN, IgG, SERUM

Test ID: ROSC

Secondary ID: 110200

EXPLANATION: The Rubeola (Measles) Antibodies, IgG, Serum test is undergoing a methodology change effective June 11, 2013.

RECOMMENDED ALTERNATE TEST: Measles (Rubeola) Antibodies IgG, Serum

TEST ID: ROPG

USEFUL FOR: Determination of immune status of individuals to the measles virus. Documentation of previous infection with measles virus in an individual with no previous record of immunization to the measles virus

NOTE: This test will contain a result code to report the Measles IgG antibody index value.

METHODOLOGY: Multiplex Flow Immunoassay

REFERENCE VALUES:

Vaccinated: Positive

Unvaccinated: Negative

SPECIMEN REQUIREMENTS:**Container/Tube:****Preferred:** Red Top**Acceptable:** Serum gel**Specimen Volume:** 0.5 mL**Minimum Volume:** 0.4 mL**SPECIMEN STABILITY INFORMATION:**

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Frozen	14 days

CPT CODE: 86765**DAY(S) SET UP:**Monday through Friday, 3rd shift**ANALYTIC TIME:**

Same day / 1 day

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or Rita M. Baird, MML Laboratory Technologist Resource Coordinator
Telephone: 800-533-1710