

**RUBEOLA (MEASLES) ANTIBODIES SCREEN, IGG, SERUM**

Test ID: ROSC

Secondary ID: 80262

**EXPLANATION:** Effective June 11, 2013, this test will be replaced with a new test code. ROPG, Measles (Rubeola) Antibodies, Serum, will contain a result code to report the Measles antibody index value.

**RECOMMENDED ALTERNATIVE TEST:** ROPG, Measles (Rubeola) Antibodies, Serum

**USEFUL FOR:** Documentation of immune status of individuals to the measles virus.  
Documentation of previous infection with measles virus in an individual without previous record of immunization to measles virus.

**METHODOLOGY:** Multiplex Flow Immunoassay (MFI)

**REFERENCE VALUES:** Vaccinated: Positive  
Unvaccinated: Negative

**SPECIMEN REQUIREMENTS:** Container/Tube: Preferred: Red top  
Acceptable: Serum gel  
Specimen Volume: 0.5 mL

**SPECIMEN STABILITY INFORMATION:**

Specimen Type	Temperature	Time
Serum	Refrigerate (Preferred)	14 days
	Frozen	14 days

**CAUTIONS:** Serum samples drawn early during acute phase of infection may be negative for IgG-class antibodies to this virus.

**CPT CODE:** 86765

**DAY(S) SET UP:** Monday through Saturday      **ANALYTIC TIME:** 1 day

**NOTE:** The following referral test code(s) will become obsolete.

Test Name	Test ID	Referral Lab Code	Referral Lab
Measles (Rubeola) IgG Ab (Immunity Screen), Serum	ZW 83 90522	964	Focus Diagnostics, Inc.

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or  
Kim J. Baker, MML Laboratory Technologist Resource Coordinator  
Telephone: 800-533-1710