

---

**Reporting Title:** Cyclosporine, B  
**Performing Location:** Rochester

**Specimen Requirements:**

Container/Tube: Lavender top (EDTA)

Specimen Volume: 3 mL

Collection Instructions:

1. Draw specimen immediately before a scheduled dose.
2. Do not centrifuge.
3. Send specimen in original tube.

Additional Information:

1. Therapeutic range applies to trough specimens drawn immediately prior to a.m. dose.
2. For specimens sent ambient, sirolimus analysis cannot be performed; if both are ordered, only cyclosporine result will be released. If sirolimus analysis is desired, specimen must be sent refrigerated or frozen.

Specimen Type	Temperature	Time
Whole Blood EDTA	Refrigerated (preferred)	14 days
	Ambient	14 days
	Frozen	14 days

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
35143	Cyclosporine, B	Numeric	ng/mL	In Process

**CPT Code:** 1 × 80158**Reference Values:**

100-400 ng/mL