

Hunter Syndrome, Known Mutation

Reporting Title: Hunter Syndrome, Known Mutation

Performing Location: Rochester

Specimen Requirements:

This test can only be performed if a mutation has previously been identified in a family member of this individual.

Forms

- 1. Molecular Genetics-Biochemical Disorders Patient Information Sheet (Supply T527) in Special Instructions
- 2. New York Clients-Informed consent is required. Please document on the request form or electronic order that a copy is on file. An Informed Consent for Genetic Testing (Supply T576) is available in Special Instructions.

Specimen must arrive within 96 hours of collection.

Submit only 1 of the following specimens:

Preferred:

Specimen Type: Whole blood

Container/Tube:

Preferred: Lavender top (EDTA) or yellow top (ACD)

Acceptable: Any anticoagulant Specimen Volume: 3 mL Collection Instructions:

Invert several times to mix blood.
 Send specimen in original tube.

Specimen Stability Information: Ambient (preferred)/Refrigerated

Specimen Type: Cultured fibroblasts Container/Tube: T-75 or T-25 flask

Specimen Volume: 1 full T-75 or 2 full T-25 flasks

Specimen Stability Information: Ambient (preferred)/Refrigerated <24 hours

Due to the complexity of prenatal testing, consultation with the laboratory is required for all prenatal testing. Prenatal specimens can be sent Monday through Thursday and must be received by 5 p.m. CST on Friday in order to be processed appropriately. All prenatal specimens must be accompanied by a maternal blood specimen. Order MCC/88636 Maternal Cell Contamination, Molecular Analysis on the maternal specimen.

Specimen Type: Amniotic fluid

Container/Tube: Amniotic fluid container

Specimen Volume: 20 mL

Specimen Stability Information: Refrigerated (preferred)/Ambient

Specimen Type: Chorionic villi

Container/Tube: 15-mL tube containing 15 mL of transport media

Specimen Volume: 20 mg

Specimen Stability Information: Refrigerated

Acceptable:

Specimen Type: Confluent cultured cells

Container/Tube: T-25 flask



Hunter Syndrome, Known Mutation

Specimen Volume: 2 flasks

Collection Instructions: Submit confluent cultured cells from another laboratory.

Specimen Stability Information: Ambient (preferred)/Refrigerated

Specimen Type: Skin biopsy

Container/Tube: Sterile container with any standard cell culture media (eg, minimal essential media, RPMI 1640). The solution should be supplemented with 1% penicillin and streptomycin. Tubes can be supplied upon request (Eagle's

minimum essential medium with 1% penicillin and streptomycin [Supply T115]).

Specimen Volume: 4-mm punch

Specimen Stability Information: Refrigerated (preferred)/Ambient

Acceptable:

Specimen Type: Blood spot

Container/Tube:

Preferred: Collection card (Whatman Protein Saver 903 Paper)

Acceptable: Ahlstrom 226 filter paper, or Supplemental Newborn Screening Card (Supply T493)

Specimen Volume: 2 to 5 Blood Spots on collection card (Whatman Protein Saver 903 Paper; Ahlstrom 226 filter paper;

or Supplemental Newborn Screening Card, Supply T493)

Collection Instructions:

- 1. An alternative blood collection option for a patient >1 year of age is finger stick.
- 2. Let blood dry on the filter paper at ambient temperature in a horizontal position for 3 hours.
- 3. Do not expose specimen to heat or direct sunlight.
- 4. Do not stack wet specimens.
- 5. Keep specimen dry

Specimen Stability Information: Ambient (preferred)/Refrigerated

Specimen Type	Temperature	Time
Varies	Varies	

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
FBC	CG030	Source: • Autopsy • Chorionic Villi • Products of Conception or Stillbirth	Answer List	No
FBC	CG033	Reason for Referral	Plain Text	No
AFC	CG132	Reason For Referral	Plain Text	No



Hunter Syndrome, Known Mutation

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
33921	Reason for Referral	Alphanumeric		In Process
33922	Result	Alphanumeric		In Process
33923	Interpretation	Alphanumeric		In Process
33924	Method	Alphanumeric		In Process
33925	Amendment	Alphanumeric		In Process
33926	Specimen	Alphanumeric		In Process
33927	Source	Alphanumeric		In Process
33928	Reviewed By	Alphanumeric		In Process
33929	Release Date	Alphanumeric		N/A

CPT Code: 1 × 81403

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
FBC	Fibroblast Culture for Genetic Test		Profile	No	Yes
AFC	Amniotic Fluid Culture/Genetic Test		Profile	No	Yes
MCC	Maternal Cell Contamination, B	1	81265	No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
FBC	16167	Specimen	Alphanumeric		31208-2
FBC	16337	Specimen ID	Alphanumeric		N/A
FBC	CG030	Source	Alphanumeric		N/A
FBC	16168	Order Date	Alphanumeric		N/A
FBC	CG033	Reason For Referral	Alphanumeric		42349-1
FBC	16172	Method	Alphanumeric		In Process
FBC	16175	Interpretation	Alphanumeric		69965-2



Hunter Syndrome, Known Mutation

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
FBC	16176	Amendment	Alphanumeric		In Process
FBC	16177	Consultant	Alphanumeric		N/A
FBC	16178	Report Date	Alphanumeric		N/A
AFC	16179	Specimen	Alphanumeric		31208-2
AFC	16338	Specimen ID	Alphanumeric		N/A
AFC	16793	Source	Alphanumeric		N/A
AFC	16180	Order Date	Alphanumeric		N/A
AFC	CG132	Reason For Referral	Alphanumeric		42349-1
AFC	16184	Method	Alphanumeric		In Process
AFC	16187	Interpretation	Alphanumeric		69965-2
AFC	16188	Amendment	Alphanumeric		In Process
AFC	16189	Consultant	Alphanumeric		N/A
AFC	16190	Report Date	Alphanumeric		N/A
мсс	23332	Specimen	Alphanumeric		31208-2
мсс	23333	Specimen ID	Alphanumeric		N/A
мсс	23334	Source	Alphanumeric		N/A
мсс	23335	Order Date	Alphanumeric		N/A
мсс	23338	Result	Alphanumeric		40704-9
мсс	23337	Amendment	Alphanumeric		N/A
мсс	23341	Reviewed By:	Alphanumeric		In Process
мсс	23339	Release Date	Alphanumeric		N/A

Reference Values:

An interpretive report will be provided.