

# **Test Definition: FHEAB**

Heroin Analysis, Whole Blood

Reporting Title: Heroin Analysis, Blood

Performing Location: Medtox Interface

## **Specimen Requirements:**

Draw blood in a grey-top NaF/oxalate tube(s). Send 7 mL of NaF/oxalate whole blood refrigerated in a plastic vial.

Specimen Type	Temperature	Time	
Whole Blood NaFl-KOx	Refrigerated (preferred)	14 days	
	Frozen	365 days	

## Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
FHEAB	Z3442	Specimen Type	Plain Text	Yes

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z3441	Heroin	Alphanumeric		
Z3442	Specimen Type	Alphanumeric		

#### **CPT Code Information:**

80101 - single drug

82542 - Quantitation (if appropriate)

## **Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
FMOR6	6-Acetylmorphine Confirm			No	No
FMORC	Morphine, unconjugated			No	No



# **Test Definition: FHEAB**

Heroin Analysis, Whole Blood

### **Result Codes for Reflex Tests:**

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
FMOR6	Z3443	6-Acetylmorphine	Alphanumeric		
FMORC	Z3444	Morphine, unconjugated	Alphanumeric		In Process

#### **Reference Values:**

This specimen was screened by Immunoassay. Any positive result is confirmed by Gas Chromatography with Mass Spectrometry (GC/MS). The following threshold concentrations are used for this analysis:

Drug Screening Threshold Confirmation Threshold

Heroin 10 ng/mL

Morphine 10 ng/mL 6-Acetylmorphine 6 ng/mL

Test Performed by: Medtox Laboratories, Inc.

402 W. County Road D St. Paul, MN 55112