

# **Test Definition: FDXAP**

Dexedrine (Dextroamphetamine)

**Reporting Title:** Dextroamphetamine **Performing Location:** Medtox Interface

### **Specimen Requirements:**

Submit only 1 of the following specimens:

#### Plasma

Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 5 mL sodium heparin plasma refrigerated in a plastic vial.

#### Serum

Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	
Varies	Refrigerated (preferred)	14 days	
	Frozen	180 days	

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z3319	Dextroamphetamine	Alphanumeric		In Process

## **CPT Code Information:**

82145

### Reference Values:

Reference Range: 10 - 100 ng/mL

Test Performed by: Medtox Laboratories, Inc.

402 W. County Road D St. Paul, MN 55112