

**REPORTING/FILE DEFINITION
CHANGE****NOTIFICATION DATE:** July 31, 2013**EFFECTIVE DATE:** September 4, 2013**PARANEOPLASTIC AUTOANTIBODY EVALUATION, SPINAL FLUID****Test ID: PNEOE**

NOTE: Effective with this change, all orders placed for this test must be submitted in the alpha-numeric character Test ID format instead of the numeric Secondary ID. Result codes may also be subject to change. Please review Test Set-Up and Conversion Mapping information at: <http://www.mayomedicallaboratories.com/test-notifications/index.html>

EXPLANATION: NMOC will be replaced as a possible reflex by NMOCC - NMO/AQP4-IgG CBA, CSF. The CBA (cell-binding) assay has a higher sensitivity than the currently offered indirect immunofluorescence assays.

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or
Steven Monson, MML Laboratory Technologist Resource Coordinator
Telephone: 800-533-1710