



MAYO CLINIC
Mayo Medical Laboratories

1-800-533-1710

**TEST UPDATE/FILE DEFINITION
CHANGE**

REFERRAL

NOTIFICATION DATE: November 15, 2013

EFFECTIVE DATE: December 16, 2013

URANIUM, URINE
Test ID: FURA
Secondary ID: 90316

EXPLANATION OF CHANGE: Test ID FURA referred to NMS Labs, will reflect the following update effective December 16, 2013.

CURRENT SPECIMEN STABILITY:

Specimen Type	Temperature	Time
Urine	Refrigerate (preferred)	28 days
	Ambient	28 days
	Frozen	28 days

NEW SPECIMEN STABILITY:

Specimen Type	Temperature	Time
Urine	Refrigerate (preferred)	14 days
	Ambient	9 days
	Frozen	30 days

QUESTIONS: Contact Mary Erath, MML Referrals Supervisor
Telephone: 800-533-1710