



MAYO CLINIC
 Mayo Medical Laboratories

1-800-533-1710

**TEST UPDATE/FILE DEFINITION
 CHANGE**

REFERRAL

NOTIFICATION DATE: November 21, 2013

EFFECTIVE DATE: Immediately

RECOMBX MATa ANTIBODY TEST

Test ID: FRECM

Secondary ID: 91892

EXPLANATION OF CHANGE: Test ID FRECM referred to Athena Diagnostics, will reflect the following update effective immediately.

CURRENT SPECIMEN REQUIRED:

Collect red-top no additive blood, separate and submit 2 mL serum in screw cap plastic vial, shipped ambient.

NEW SPECIMEN REQUIRED:

Collect red-top no additive blood, separate and submit 2 mL serum in screw cap plastic vial, shipped refrigerate.

CURRENT SPECIMEN STABILITY:

Specimen Type	Temperature	Time
Serum Red	Ambient (preferred)	
	Refrigerate	

NEW SPECIMEN STABILITY:

Specimen Type	Temperature	Time
Serum Red	Refrigerate (preferred)	14 days
	Ambient	72 hours
	Frozen	

QUESTIONS: Contact Mary Erath, MML Referrals Supervisor
 Telephone: 800-533-1710