**Test Definition: FMT21** 

MaterniT21 Plus

Reporting Title: MaterniT21 Plus

Performing Location: Sequenom Center for Molecular Medicine LLC

### **Specimen Requirements:**

Specimen should be collected between 10 and 27 weeks' gestation.

Sequenom collection kit (MML Supply T715) required:

Draw 2 tubes of blood, 10 mL each in special Streck Black/Tan top tube kit (MML supply number T715). Ship ambient

Specimen Type	Temperature	Time
WB Streck	Ambient	72 hours

## Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
FMT21	Z3731	Maternal Height (inches)	Plain Text	No
FMT21	Z3732	Maternal Weight (pounds)	Plain Text	No
FMT21	Z3733	Gestational age (weeks)	Plain Text	No
FMT21	Z3734	Gestational age (days)	Plain Text	No
FMT21	Z3735	Method determining gestational age	Plain Text	No
FMT21	Z3736	Number of fetuses	Plain Text	No
FMT21	Z3737	Increased risk due to	Plain Text	No

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z3738	Result	Alphanumeric		In Process

**CPT Code:** 1 × 81479



# **Test Definition: FMT21**

MaterniT21 Plus

#### **Reference Values:**

A final report will be attached in MayoAccess

Test Performed by: Sequenom Center for Molecular Medicine LLC 3595 John Hopkins Court San Diego, CA 92121