

TEST DEFINITION

2/16/2009

Code name

83129 Methadone Confirmation, U

**** End of Name to order code Report ****

MML		MML Test setup information				Perform
ORDER	EFF	TC	TITLE	Checking Normals	Print normals (# coded)	Site *
CODE	DATE					
-----	-----	-----	-----	-----	-----	-----
83129	9/26/2008		Methadone Confirmation, U			MCR
			Transport temp : Refrig\Frozen OK\Ambient <72 hours OK			
		21105	Immunoassay Screen			
			NO SEX			
			All Ages :		; Cutoff: 300	
			MALE			
			All Ages :		; Cutoff: 300	
			FEMALE			
			All Ages :		; Cutoff: 300	
			- - - - -			
		56028	EDDP			
			Possible result values include : Negative			
			Units: ng/mL			
			NO SEX			
			All Ages :		; Cutoff: 100	
			MALE			
			All Ages :		; Cutoff: 100	
			FEMALE			
			All Ages :		; Cutoff: 100	
			- - - - -			
		83129	Methadone-by GC/MS			
			Expect results of form : Should be a 1, 2, 3, 4, or 5 digit whole number			
			Possible result values include : Negative			
			Units: ng/mL			
			NO SEX			
			All Ages :		; Cutoff: 100	
			MALE			
			All Ages :		; Cutoff: 100	

(continued next page)

MML MML Test setup information

ORDER CODE -----	EFF DATE -----	TC -----	TITLE -----	Checking Normals -----	Print normals (# coded) -----	Perform Site * -----
83129	(continued)					
		83129	Methadone-by GC/MS (continued)			
			FEMALE			
			All Ages :		; Cutoff: 100	
			- - - - -			
		21107	Interpretation			
			- - - - -			
		21110	Chain of Custody			
			Possible result values include : ON FILE			
			- - - - -			

*Performing Site Legend

```

=====
==
MCR      Mayo Clinic Dpt of Lab Med & Pathology      LAB DIRECTOR:  Franklin R. Cockerill, III, M.D.
          200 First Street SW
          Rochester, MN 55905
=====
==

```

*** End of Report ***

MML Messages used as normals

```

CODE      TEXT
-----

```

Total of 0 normals codes

*** End of Report ***

To obtain test set-up information contact MML Client Support or log on to:
<http://www.mayomedicallaboratories.com/customer-service/set-up/index.html>